BENEFITS

Delicately delivers treatments for improving skin quality and skin health

- Pain-free delivery
- Direct dermal application
- Easy, quick, and simple procedure

Painless treatment of sensitive areas

- Perioral (around the mouth)
- Periocural (around the eyes)
- Forehead
- Neck
- · Back of hands
- Décolleté

0.18 mm 0.13 mm Hair Microchannel

Professional Model channel length: 600μ, 1000μ Consumer Model channel length: 250μ

Acne Scars

Before

48 hours post treatment with Hyaluronic Acid







Periocular

Before

7 days post treatment (Micro Botox®)







Histology test conducted by Dr. Daniel Del Vecchio and Dr. Louis Bucky during Lipoplasty on live human skin.

A Q U A G O L D[®]

fine touch™

Exclusively available in your physician's office for provider-quided treatments

Experience clinical luxury following intradermal treatments, including:

- Hyaluronic acid fillers
- Neuromodulators
- Vitamin complex and mesotherapy
- Growth factors and PRP

What our doctors have to say about AQUAGOLD® fine touch™

"Unlike classic roller techniques that can cause microtears in the skin, this device won't cause bleeding because it applies the microchannel perpendicular to the skin."

> Joel Schlessinger, MD Board-Certified Dermatologist (Omaha, NE)

"What I like about the Aquagold is that you have the versatility to put whatever content you would like in the chamber to increase penetration into the skin."

> Nathan Trookman, MD Board-Certified Dermatologist (Colorado Springs, CO)

"This is a very efficient device where a reservoir can be used to introduce substances like HA, botulinum toxin and other substances that may improve skin quality in combination."

> Steven Fagien, MD Board-Certified Oculoplastic Surgeon (Boca Raton, FL)



Painlessly, easily, and effectively delivers treatments with 99.9% pure gold-plated surgical grade stainless steel microchannels, each thinner than a human hair.

"Superficial microinjections of diluted HA using a new microchannel device appears to provide significant visual improvement in moderate to severe acne scarring and skin quality with no adverse effects."

Excerpt from Study presented by Erin Gilbert, MD. PhD at South Beach Symposium 2016

AS FEATURED IN:



For international inquiries call 702-874-8544









Copyright © 2014-2016 Aquavit Pharmaceuticals, Inc. Science Personalized, Healthcare Modernized.

> All registered trademarks are the property of their respective owners. AQTMARCOL892592004



fine touch™

DEFINE YOUR CLINICAL LUXURY™ WITH NEW MICRO-CHANNELING TECHNOLOGY



www.aquagoldfinetouch.com 1-800-55-fine touch (1-800-553-4638)

info@aquagoldfinetouch.com 1-800-55-fine touch (1-800-553-4638)

Ask about our Loyalty Program: • Heart of Gold for MDs • Gold Rush for patients